


Energy nonlinearity correction of a high resolution metallic magnetic calorimeter



M.L. Zahir, A. Kaur, M. Loidl, M. Rodrigues

Université Paris-Saclay, CEA, List, Laboratoire National Henri Becquerel (LNE-LNHB), F-91120, Palaiseau, France

Precise measurements of X-ray and Y-ray energies are important for the study of atomic and nuclear structure. They contribute to the advancement of our fundamental knowledge in physics. Cryogenic detectors offer excellent energy resolution and are ideal for precise energy measurements, but they have intrinsic and extrinsic nonlinearities. It is crucial for cryogenic detector to correct their nonlinearity in order to accurately interpret the measured energies. The main nonlinearity is intrinsic to the detector. Other contributions to nonlinearity can arise from the electronics and the signal processing. Therefore, it is important to assess the impact of these contributions for the subsequent nonlinearity



Conclusion

59.318847(50)

Tungsten

The detector nonlinearity remains predominant. We observe a small impact on the linearity due to signal processing.

0.025

For now, the gain does not have any impact on the nonlinearity witch mean that we do not observe nonlinearity from the electronics. However if in the future we are able to correct the intrinsic nonlinearity of the detector, it is possible that the electronics becomes an source on the nonlinearity. The theoretical calculation of the nonlinearity is consistent with the experimental data using the optimal method. This shows that the calculations are reliable enough to be used for correcting the nonlinearity in subsequent measurement. However, it is still important to confirm the calculations reliability using additional data.

The full template and the area methods have a better

linearity than the optimal filter and the short template