

# Immobilized nitrilase in continuous bioreactors: packed bed reactor vs. milli-structured pillar reactor

Maïté Michaud, Carine Vergne-Vaxelaire, Véronique de Berardinis, Anne Zaparucha, Guillaume Nonglaton, Zoe Anxionnaz-Minvielle

#### ▶ To cite this version:

Maïté Michaud, Carine Vergne-Vaxelaire, Véronique de Berardinis, Anne Zaparucha, Guillaume Nonglaton, et al.. Immobilized nitrilase in continuous bioreactors: packed bed reactor vs. millistructured pillar reactor. ECCEECAB 2023 - 14th European Congress of Chemical Engineering and 7th European Congress of Applied Biotechnology, Sep 2023, Berlin, Germany. cea-04431029

# HAL Id: cea-04431029 https://cea.hal.science/cea-04431029

Submitted on 1 Feb 2024

**HAL** is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers.

L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.

# Immobilized nitrilase in continuous bioreactors: packed bed reactor vs. milli-structured pillar reactor

Maïté MICHAUD<sup>a</sup>, Carine VERGNE-VAXELAIRE<sup>c</sup>, Véronique DE BERARDINIS<sup>c</sup>,
Anne ZAPARUCHA<sup>c</sup>, Guillaume NONGLATON<sup>b</sup>, Zoé ANXIONNAZ-MINVIELLE<sup>a</sup>

<sup>a</sup> Univ. Grenoble Alpes, CEA, LITEN, DTCH, F-38000, Grenoble, France

<sup>b</sup> Univ. Grenoble Alpes, CEA, LETI, DTBS, F-38000, Grenoble, France

<sup>c</sup> Génomique Métabolique, Genoscope, Institut François Jacob, CEA, CNRS, Univ.

Evry, Univ. Paris-Saclay, 91057 Evry, France

## **Biocatalysis flow context**

In order to address the sustainability and energy efficiency issues, biocatalytic synthesis routes have increasingly gained more popularity over the years. This green tool has many benefits to offer such as soft operating conditions (i.e low temperature and pressure), high selectivity, and aqueous medium reducing the demand in petrochemical solvents [1]. To make industrial investments more cost-effective, immobilization strategies are actively studied by the scientific community to propose efficient and stable immobilization techniques. Another profitable aspect of immobilized enzymes is the possible use of continuous processes instead of traditional batch-wise production. However, enzyme immobilization is far from being a mature discipline. Thus, process intensification from batch to continuous production is one of the main point to address for successful deployment of flow biocatalysis at an industrial level [2].

#### Scope of the study and objectives

A promising primary approach have been published with continuous production of nicotinic acid over 30 days (i.e. ≈ 1400 cycles) using a traditional Packed-Bed-Reactor "PBR" with PMMA beads <sup>[3]</sup>. Considering the ending conversion of 50%, the total turnover "TN" (i.e mass of product per mass of immobilized enzyme) was 8-fold higher with PBR compared to cyclic use of the immobilized enzymes in batch conditions. Even higher performances were obtained changing from a commercial to a thermostable nitrilase from an in-house collection with a fully stable conversion of 100%. However, PBR are known for excessive back-pressure which hampers the scale-up feasibility. Surface coated-reactors may remove this barrier but reported datas in the literature are focused on micro-equipments. For an industrial scale application, both scale-up (i.e. increasing the characteristic dimensions of the channel) and numbering-up (i.e.

increasing the number of reactor) must be considered <sup>[4]</sup>. For this study, we immobilized the commercial nitrilase in a milli-scale pillar reactor (i.e. channel length of 1mm). The reactor walls were coated with SiO<sub>2</sub> beads for higher loading capacity.

### **Biocatalytic performances**

The experimental results (Figure 1) show a lower stability in the milli-reactor than in the PBR. This may be due to the enzyme immobilization protocol differences ( $SiO_2$  vs. PMMA beads). This aspect will be studied later for the reactor optimization. Compared to the batch reactor using PMMA beads, with similar ending conversion at 50%, the productivity capacity is greatly improved with the continuous mode: 14 µmol/h instead of 1 µmol/h. The TN is also 2-fold higher with milli-structured reactor compared to the batch mode.




Figure 1: (Left) 3D CAO-view of the milli-pillar reactor; (Right) Biocatalytic conversion of the PBR and millireactor - Both reactors have an internal volume of 0.5 ml and the liquid flowrate is 0.03 ml/min

## **Conclusion and perspectives**

Continuous biocatalytic synthesis with a milli-pillar reactor gives greater TN and productivity compared to a batch mode production. The process scalability potential is enhanced changing from PBR to a milli-structured reactor by a decrease in the pumping power (i.e theoretical pressure drop at least 50% higher for PBR than the one measured with millireactor). Multiple scale-up stages are planned for a final aiming flow of 1L/h by combining CFD modeling and experimental tests. Also, the reactor structure influence will be considered with comparison with a zigzag meandering milli-channel.

- [1] G. Rossino et al., Chirality, 34 (11), 2022, 1403-1418.
- [2] J.M. Bolivar et al., Chem. Soc. Rev., 51, 2022, 6251-6290.
- [3] C. Teepakorn et al., *Biotechnol.J.*, 2021;16:2100010.
- [4] S.G. Newman and K.F. Jensen, *Green Chem.*, 15, 2013, 1456-1472.